EEG Processing Report

Run ID: 01K4AHHYHQR91EZRXSCD27PFB3 - SUCCESS

Generated: 2025-09-04 10:22:39

Import Information Processing Details 128_Rest_EyesOpen_D1004.raw File Filter 1-80 Hz Duration 303.4 sec Notch 60, 120 Hz Sample Rate 1000.0 Hz Resample Rate 250 Hz Channels 128 Trim Duration 4 sec Reference average

Task Parameters

Parameter	Value
move_flagged_files	False
resample_step.enabled	True
resample_step.value	250
filtering.enabled	True
filtering.value.l_freq	1
filtering.value.h_freq	80
filtering.value.notch_freqs	[60, 120]
filtering.value.notch_widths	5
drop_outerlayer.enabled	False
drop_outerlayer.value	0
eog_step.enabled	False
eog_step.value	[1, 32, 8, 14, 17, 21, 25, 125, 126, 127, 128]
trim_step.enabled	True
trim_step.value	4
crop_step.enabled	False
crop_step.value.start	0
crop_step.value.end	60
reference_step.enabled	True
reference_step.value	average
montage.enabled	True
montage.value	GSN-HydroCel-129
ICA.enabled	True
ICA.value.method	infomax
ICA.value.n_components	None
ICA.value.fit_params.extended	True
component_rejection.enabled	True
component_rejection.method	icvision
component_rejection.value.ic_flags_to_reject	['muscle', 'heart', 'eog', 'ch_noise', 'line_noise']

component_rejection.value.ic_rejection_threshold	0.3
epoch_settings.enabled	True
epoch_settings.value.tmin	-1
epoch_settings.value.tmax	1
epoch_settings.event_id	None
epoch_settings.remove_baseline.enabled	False
epoch_settings.remove_baseline.window	[None, 0]
epoch_settings.threshold_rejection.enabled	False
epoch_settings.threshold_rejection.volt_threshold.eeg	0.000125

Processing Steps

Processing Step		
Create Bids Path		
Resample Data		
Filter Data		
Trim Edges		
Clean Bad Channels		
Rereference Data		
Annotate Noisy Epochs		
Annotate Uncorrelated Epochs		
Detect Dense Oscillatory Artifacts		
Run Ica		
Apply Ica Component Rejection		
Create Regular Epochs		
Detect Outlier Epochs		
Gfp Clean Epochs		
Psd Topo Figure		

Bad Channels

Source	Bad Channels
Clean Bad Channels	E3, E122, E107, E8, E14, E56, E17, E127, E2
Uncorrelated Channels	E2, E3, E17, E56, E107, E122
Deviation Channels	E8, E14, E17
Ransac Channels	E127
Total Removed	9 / 128 (7.0%)

Results Summary

Metric	Value	
Processing State	completed_files	
Epochs Retained	90 / 147 (61.2%)	
Duration Retained	180.0s / 294.0s (61.2%)	
Removed ICA Components	0, 1, 2, 4, 7, 8, 12, 13, 17, 18, 20, 28, 30, 31, 32, 35, 38, 45, 50, 64,	70, 98, 11

Quality Control Flags

Flagged Reason

WARNING: 130 potential reference artifacts detected

Output Files

File Name 128_Rest_EyesOpen_D1004_psd_topo_figure.png 128_Rest_EyesOpen_D1004-ica.fif 128_Rest_EyesOpen_D1004_raw_vs_cleaned_overlay.png sub-451559_task-restingstatetutorialica_components_all.pdf

Run ID: 01K4AHHYHQR91EZRXSCD27PFB3 | Task: RestingStateTutorial | Timestamp: N/A